

Measurements for Selecting Wunzi Size (Metric)

Torso Length (in CENTIMETERS) (see Fig. 1) - Measure down from C-7 (the big bone at the bottom of the neck) over the diaper and up to the child's navel. Be careful to wrap the tape gently under the diaper.

Belly Girth - Measure the infant's belly circumference at the navel (in CENTIMETERS).

NOTE: AGE AND WEIGHT ARE GUIDELINES ONLY - DO NOT ORDER WITHOUT THE TORSO AND GIRTH MEASUREMENTS!

Model #	Model Name	Age	Weight	Torso Length	Belly Girth
GWZ25	Wunzi XXSmall	0 – 3 mo.	3.2 - 5.4 kg	40 cm	40 cm
GWZ50	Wunzi XSmall	3 – 6 mo.	5.4 - 7.7 kg	50 cm	44 cm
GWZ75	Wunzi Small	6 – 12 mo.	7.7 - 8.2 kg	55 cm	46 cm
GWZLKS	Limb Kit, Small	0 -12 mo.	Up to 8.2 kg		
GWZ100	Wunzi Medium	12 – 18 mo.	8.2 - 9.5 kg	60 cm	50 cm
GWZ125	Wunzi Large	18 – 24 mo.	9.5 - 11.8 kg	68 cm	54 cm
GWZLKL	Limb Kit, Large	12 – 24 mo.	Up to 11.8 kg		



(Fig. 1)Torso Length